


Rugged Construction! Corrosion Resistant!



MASTER-POX

Convenient
Storage
Case (included)



NEW Durable Material!

Ideal for measuring salty and acidic samples

MASTER series

MADE IN JAPAN

MASTER series Refractometer
(Automatic Temperature Compensation & Water Resistant)

MASTER-POX Cat.No. 2391



**Water
Resistant (IP65)**


**Automatic
Temperature
Compensation**

scale

Scale range	Brix 0.0~33.0%
Minimum Scale	Brix 0.2%
Accuracy	Brix $\pm 0.2\%$ (10~30°C)
Repeatability	Brix $\pm 0.1\%$
International Protection Class	IP65 (except eye piece) Dust-tight and Protected against water jets.
Size & Weight	3.3 × 3.3 × 20.4cm, 100g

MASTER series Refractometer

MASTER-PM Cat.No. 2393



scale

Scale range	Brix 0.0~33.0%
Minimum Scale	Brix 0.2%
Size & Weight	3.3 × 3.3 × 20.4cm, 100g



All ATAGO refractometers are designed and manufactured in Japan.

HACCP | G.M.P. | G.L.P.

ATAGO products are acknowledged by HACCP, G.M.P. and G.L.P.

ATAGO CO., LTD.

Headquarters: 32-10, Honcho, Itabashi-ku, Tokyo 173-0001 Japan
TEL: 81-3-3964-6156 FAX: 81-3-3964-6137
overseas@atago.net http://www.atago.net/

ATAGO U.S.A., Inc.

12011 NE First Street, Bldg. C, Suite 110, Bellevue, WA 98005 U.S.A.
TEL: 1-425-637-2107 FAX: 1-425-637-2110
customerservice@atago-usa.com http://www.atago.net/

ATAGO INDIA Instruments Pvt. Ltd.

A-501, Mangal Aarambh Business Centre, Near Kora Kendra, Off. S.V. Road, Borivali (West), Mumbai-400 092 India
TEL: 91-22-2833-8038 / 8076 FAX: 91-22-2899-8164
customerservice@atago-india.com



※Specifications and appearance are subject to change without notice.

0807K

HAND-HELD REFRACTOMETERS