

DIGITAL BUTYRO  
REFRACTOMETER

# PR-BUTYRO

Cat.No.3454

New ATAGO Digital Model PR-BUTYRO perfect for Butter, Ghee & Edible Oil industry for checking purity / adulteration.

## Features

- Calibration with Standard liquid
- Large, easy to read display
- External Light Interference (ELI)\*
- Display of both Butyro & RI scale
- Temperature Compensated Scale at 40°C
- Automatic & Accurate (Accuracy : Butyro ±0.5)



Dual scale of Butyro & Refractive Index !



## Specifications

Measurement Range (Automatic Temperature Compensation)	Butyro 30.0~90.0	RI 1.4450~1.4850 (Converted at 40°C)
Minimum Indication	Butyro 0.1	RI 0.0001
Measurement Accuracy	Butyro ±0.5 (at 40°C)	RI ±0.0003 (at 40°C)
Measurement Temperature	10~50°C (Automatic Temperature Compensation)	
Ambient Temperature	10~40°C	
Sample Volume	0.1ml or more	
Measurement Time	3 seconds (approximately)	
Power Supply	006P dry battery (9V)	
International Protection Class	IP 64 Water Resistant	
Dimensions and Weight	17(W) × 9(D) × 4(H)cm, 300g (Main unit only)	

### \*External Light Interference (ELI)

When intense light penetrates the prism of a digital refractometer, the light waves may interfere with the sensor resulting in inaccurate readings.

To ensure accurate measurement results each time, the PR-BUTYRO is programmed with the ELI function to display the warning message "nnn" when intense direct light is detected. Forming a habit of shading the sample stage with your hand and re-pressing the START/OFF Switch (when the warning message from the ELI function is displayed) will ensure accurate measurement results each time.

All ATAGO refractometers are designed and manufactured in Japan.



ATAGO products are acknowledged by HACCP, G.M.P. and G.L.P.

### ATAGO CO.,LTD.

Headquarters: 32-10, Honcho, Itabashi-ku, Tokyo 173-0001 Japan  
TEL: 81-3-3964-6156 FAX: 81-3-3964-6137  
overseas@atago.net http://www.atago.net/

### ATAGO U.S.A.,Inc.

12011 NE First Street, Bldg. C, Suite 110, Bellevue, WA 98005 U.S.A.  
TEL: 1-425-637-2107 FAX: 1-425-637-2110  
customerservice@atago-usa.com http://www.atago.net/

### ATAGO INDIA Instruments Pvt. Ltd.

A-501, Mangal Aarambh Business Centre, Near Kora Kendra, Off. S.V. Road, Borivali (West), Mumbai-400 092 India  
TEL: 91-22-2833-8038 / 8076 FAX: 91-22-2899-8164  
customerservice@atago-india.com



※Specifications and appearance are subject to change without notice.

