

DIGITAL HAND-HELD "POCKET" REFRACTOMETER

65th anniversary

PAL- α

Cat.No.3840

Special Anniversary model
with a wide Brix range 0 to 85%!!



Actual Size

Application

- * Outdoors in the vegetable fields and fruit orchards
- * Indoors for beverages (fruits juice, coffee, soft drink, etc) and foods (soy sauce, sauce, seasoning, soup, jam, marmalade, maple syrup, etc.)
- * At fruit and vegetable packing facilities and in the supermarket
- * At restaurants or shops
- * At junior high schools, high schools, or universities as course material

Features

- * Calibration with water only
- * Measurement in 3 seconds
- * Extremely water resistant (IP65)
- * Light & Compact, 100g
- * Ergonomically designed for easy, one-handed operation
- * Convenient storage case
- * Automatic Temperature Compensation (ATC)
- * Designed with a strap holder option



Measurement Method



Apply 2-3 drops onto the prism surface.



Press the START key.



The measurement value will be displayed in 3 seconds.

Specifications

Measurement range	: Brix 0.0 to 85.0%
Resolution	: Brix 0.1%
Measurement accuracy	: Brix $\pm 0.2\%$
Repeatability	: Brix $\pm 0.1\%$
Ambient Temperature	: 10 to 40°C
Measurement Temperature	: 10 to 75°C (Automatic Temperature Compensation)
Sample volume	: 0.3ml or more
Measurement time	: 3 seconds
Power supply	: AAA alkaline battery $\times 2$
Battery life	: Approx. 11,000 measurements (Using new alkaline batteries)
International Protection Class	: IP65 Water resistant
Dimensions and weight	: 55(W) \times 31(D) \times 109(H)mm, 100g

All ATAGO refractometers are designed and manufactured in Japan.

HACCP GMP GLP

ATAGO products comply with HACCP, GMP, and GLP system standards.

ATAGO CO., LTD.

Headquarters: The Front Tower Shiba Koen, 23rd Floor
2-6-3 Shiba-koen, Minato-ku, Tokyo 105-0011, Japan
TEL : 81-3-3431-1943 FAX : 81-3-3431-1945
overseas@atago.net http://www.atago.net/

ATAGO U.S.A., Inc.

TEL : 1-425-637-2107 customerservice@atago-usa.com

ATAGO INDIA Instruments Pvt. Ltd.

TEL : 91-22-2833-8038 / 8076 customerservice@atago-india.com

ATAGO (THAILAND) Co., Ltd.

TEL : 662-982-8718-9 customerservice@atago-thailand.com

ATAGO BRASIL Ltda.

TEL : 55 16 3916-6000 customerservice@atago-brasil.com

ATAGO ITALIA s.r.l.

TEL : 39 2 36557267 customerservice@atago-italia.com

ATAGO CHINA Guangzhou Co., Ltd.

TEL : 86 20-38106065 info@atago-china.com



* Specifications and appearance are subject to change without notice.