

Vegetable & Fruit Refractometer

MASTER-AGRI

For small-volume measurement

Optimized for measuring concentrations (Brix) of plant juices from leaves, stems and roots.

Only a small amount of liquid is needed for measurement.

Perfectly apt to measure fruit juices and more.



Optional

SUPER EXTRACTOR
Part No. RE-29401

MASTER-AGRI Cat.No.2462

Scale range	Brix 0.0 to 53.0%
Minimum Scale	Brix 0.5%
Dimensions & Weight	3.2×3.4×16.8cm, 130g

All ATAGO refractometers are designed and manufactured in Japan.



Headquarters: The Front Tower Shiba Koen,
23rd Floor 2-6-3 Shiba-koen, Minato-ku,
Tokyo 105-0011, Japan
TEL : 81-3-3431-1943 FAX : 81-3-3431-1945



ATAGO products comply with HACCP, GMP, and GLP system standards.

<http://www.atago.net/> overseas@atago.net

- ATAGO U.S.A., Inc.** TEL : 1-425-637-2107 customerservice@atago-usa.com
- ATAGO INDIA Instruments Pvt. Ltd.** TEL : 91-22-2833-8038 / 8076 customerservice@atago-india.com
- ATAGO (THAILAND) Co., Ltd.** TEL : 662-982-8718-9 customerservice@atago-thailand.com
- ATAGO BRASIL Ltda.** TEL : 55 16 3916-6000 customerservice@atago-brasil.com
- ATAGO ITALIA s.r.l.** TEL : 39 2 3657267 customerservice@atago-italia.com
- ATAGO CHINA Guangzhou Co., Ltd.** TEL : 86 20-38106065 info@atago-china.com

* Specifications and appearance are subject to change without notice.