

温度 (°C)	蔗糖濃度測定値 (質量分率)																	
	0	5	10	15	20	25	30	35	40	45	50	55	60	65	70	75	80	85
15	-0.29	-0.30	-0.32	-0.33	-0.34	-0.35	-0.36	-0.37	-0.37	-0.38	-0.38	-0.38	-0.38	-0.38	-0.38	-0.37	-0.37	
16	-0.24	-0.25	-0.26	-0.27	-0.28	-0.28	-0.29	-0.30	-0.30	-0.30	-0.31	-0.31	-0.31	-0.31	-0.31	-0.30	-0.30	
17	-0.18	-0.19	-0.20	-0.20	-0.21	-0.21	-0.22	-0.22	-0.23	-0.23	-0.23	-0.23	-0.23	-0.23	-0.23	-0.23	-0.22	
18	-0.12	-0.13	-0.13	-0.14	-0.14	-0.14	-0.15	-0.15	-0.15	-0.15	-0.15	-0.15	-0.15	-0.15	-0.15	-0.15	-0.15	
19	-0.06	-0.06	-0.07	-0.07	-0.07	-0.07	-0.07	-0.08	-0.08	-0.08	-0.08	-0.08	-0.08	-0.08	-0.08	-0.08	-0.07	
20	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	
21	+0.06	+0.07	+0.07	+0.07	+0.07	+0.07	+0.08	+0.08	+0.08	+0.08	+0.08	+0.08	+0.08	+0.08	+0.08	+0.08	+0.07	
22	+0.13	+0.14	+0.14	+0.14	+0.15	+0.15	+0.15	+0.15	+0.16	+0.16	+0.16	+0.16	+0.16	+0.16	+0.15	+0.15	+0.15	
23	+0.20	+0.21	+0.21	+0.22	+0.22	+0.23	+0.23	+0.23	+0.23	+0.24	+0.24	+0.24	+0.24	+0.23	+0.23	+0.23	+0.22	
24	+0.27	+0.28	+0.29	+0.29	+0.30	+0.30	+0.31	+0.31	+0.31	+0.32	+0.32	+0.32	+0.32	+0.31	+0.31	+0.31	+0.30	
25	+0.34	+0.35	+0.36	+0.37	+0.38	+0.38	+0.39	+0.39	+0.40	+0.40	+0.40	+0.40	+0.40	+0.39	+0.39	+0.38	+0.37	
26	+0.42	+0.43	+0.44	+0.45	+0.46	+0.46	+0.47	+0.47	+0.48	+0.48	+0.48	+0.48	+0.48	+0.47	+0.47	+0.46	+0.45	
27	+0.50	+0.51	+0.52	+0.53	+0.54	+0.55	+0.55	+0.56	+0.56	+0.56	+0.56	+0.56	+0.56	+0.55	+0.55	+0.54	+0.52	
28	+0.58	+0.59	+0.60	+0.61	+0.62	+0.63	+0.64	+0.64	+0.64	+0.65	+0.65	+0.64	+0.64	+0.63	+0.63	+0.62	+0.60	
29	+0.66	+0.67	+0.68	+0.70	+0.71	+0.71	+0.72	+0.73	+0.73	+0.73	+0.73	+0.73	+0.72	+0.72	+0.71	+0.70	+0.67	
30	+0.74	+0.76	+0.77	+0.78	+0.79	+0.80	+0.81	+0.81	+0.82	+0.82	+0.81	+0.81	+0.80	+0.80	+0.79	+0.78	+0.75	
31	+0.83	+0.84	+0.85	+0.87	+0.88	+0.89	+0.89	+0.90	+0.90	+0.90	+0.90	+0.89	+0.89	+0.88	+0.87	+0.86	+0.82	
32	+0.92	+0.93	+0.94	+0.96	+0.97	+0.98	+0.98	+0.99	+0.99	+0.99	+0.99	+0.98	+0.97	+0.96	+0.95	+0.93	+0.90	
33	+1.01	+1.02	+1.03	+1.05	+1.06	+1.07	+1.07	+1.08	+1.08	+1.08	+1.07	+1.07	+1.06	+1.04	+1.03	+1.01	+0.98	
34	+1.10	+1.11	+1.13	+1.14	+1.15	+1.16	+1.16	+1.17	+1.17	+1.16	+1.16	+1.15	+1.14	+1.13	+1.11	+1.09	+1.05	
35	+1.19	+1.21	+1.22	+1.23	+1.24	+1.25	+1.25	+1.26	+1.26	+1.25	+1.25	+1.24	+1.23	+1.21	+1.19	+1.17	+1.13	
36	+1.29	+1.30	+1.31	+1.33	+1.34	+1.34	+1.35	+1.35	+1.35	+1.34	+1.34	+1.33	+1.31	+1.29	+1.28	+1.25	+1.20	
37	+1.39	+1.40	+1.41	+1.42	+1.43	+1.44	+1.44	+1.44	+1.44	+1.43	+1.43	+1.41	+1.40	+1.38	+1.36	+1.33	+1.28	
38	+1.49	+1.50	+1.51	+1.52	+1.53	+1.53	+1.54	+1.54	+1.53	+1.53	+1.52	+1.52	+1.48	+1.46	+1.44	+1.42	+1.36	
39	+1.59	+1.60	+1.61	+1.62	+1.63	+1.63	+1.63	+1.63	+1.63	+1.62	+1.61	+1.59	+1.57	+1.55	+1.52	+1.50	+1.43	
40	+1.69	+1.70	+1.71	+1.72	+1.73	+1.73	+1.73	+1.73	+1.72	+1.71	+1.70	+1.68	+1.66	+1.63	+1.61	+1.58	+1.51	