

" 温度 °C"	pH 値				
	しゅう酸塩	フタル酸塩	中性りん酸塩	ほう酸塩	炭酸塩 (1)
0	1.67	4.01	6.98	9.46	10.32
5	1.67	4.01	6.95	9.39	(10.25)
10	1.67	4.00	6.92	9.33	10.18
15	1.67	4.00	6.90	9.27	(10.12)
20	1.68	4.00	6.88	9.22	(10.07)
25	1.68	4.01	6.86	9.18	10.02
30	1.69	4.01	6.85	9.14	(9.97)
35	1.69	4.02	6.84	9.10	(9.93)
38	-	-	-	-	9.91
40	1.70	4.03	6.84	9.07	-
45	1.70	4.04	6.83	9.04	-
50	1.71	4.06	6.83	9.01	-
55	1.72	4.08	6.84	8.99	-
60	1.73	4.10	6.84	8.96	-
70	1.74	4.12	6.85	8.93	-
80	1.77	4.16	6.86	8.89	-
90	1.80	4.20	6.88	8.85	-
95	1.81	4.23	6.89	8.83	-