

" 温度 °C"	pH 値										
	しゅう酸塩		フタル酸塩		中性りん酸塩		りん酸塩		ほう酸塩		炭酸塩
	第 1 種	第 2 種	第 1 種	第 2 種	第 1 種	第 2 種	第 1 種	第 2 種	第 1 種	第 2 種	炭酸塩
0	1.666	1.67	4.003	4.00	6.984	6.98	7.534	7.53	9.464	9.46	第 2 種
5	1.668	1.67	3.999	4.00	6.951	6.95	7.500	7.50	9.395	9.40	10.32
10	1.670	1.67	3.998	4.00	6.923	6.92	7.472	7.47	9.332	9.33	10.24
15	1.672	1.67	3.999	4.00	6.900	6.90	7.448	7.45	9.276	9.28	10.18
20	1.675	1.68	4.002	4.00	6.881	6.88	7.429	7.43	9.225	9.22	10.12
25*(2)	1.679*	1.68*	4.008*	4.01*	6.865*	6.86*	7.413*	7.41*	9.180*	9.18*	10.06
30	1.683	1.68	4.015	4.02	6.853	6.85	7.400	7.40	9.139	9.14	10.01*
35	1.688	1.69	4.024	4.02	6.844	6.84	7.389	7.39	9.102	9.10	9.97
38	1.691	1.69	4.030	4.03	6.840	6.84	7.384	7.38	9.081	9.08	9.92
40	1.694	1.69	4.035	4.04	6.838	6.84	7.380	7.38	9.068	9.07	-
45	1.700	1.70	4.047	4.05	6.834	6.83	7.373	7.37	9.038	9.04	9.89
50	1.707	1.71	4.060	4.06	6.833	6.83	7.367	7.37	9.011	9.01	9.86
55	1.715	1.72	4.075	4.08	6.834	6.83	-	-	8.985	8.98	9.83
60	1.723	1.72	4.091	4.09	6.836	6.84	-	-	8.962	8.96	-
70	1.743	1.74	4.126	4.13	6.845	6.84	-	-	8.921	8.92	-
80	1.766	1.77	4.164	4.16	6.859	6.86	-	-	8.885	8.88	-
90	1.792	1.79	4.205	4.20	6.877	6.88	-	-	8.850	8.85	-
95	1.806	1.81	4.227	4.23	6.886	6.89	-	-	8.833	8.83	-