

" 温度 °C"	0.1mol / l 水酸化 ナトリウム溶液	飽和水酸化 カルシウム溶液
0	13.8	13.43
5	13.6	13.21
10	13.4	13.00
15	13.2	12.81
20	13.1	12.63
25	12.9	12.45
30	12.7	12.30

" 温度 °C"	0.1mol / l 水酸化 ナトリウム溶液	飽和水酸化 カルシウム溶液
35	12.6	12.14
40	12.4	11.99
45	12.3	11.84
50	12.2	11.70
55	12.0	11.58
60	11.9	11.45